

Bilsler®



Since 1980, thousands of references ...

SURGICAL SUCTION UNIT



MODEL
BA-3



MODEL
BA-6

2 year guarantee

Bilsler Brand Surgical Suctions, all transactions that may require surgical interventions and aspirator; operating room, intensive care, emergency units and is designed to be used in all clinics.

MODEL	Capacity	Jar	Jar Type	Pump	Vacuum	Type	Power
BA-3	3 lt	Single	Glass with safety float	Oil free	60 lt/min 700 mmHg with vacuum adjustment control switch	Mobile 4 antistatic wheels, 2 with brakes	220-230 VAC 50/60 Hz
BA-6	2x3 lt	Double with left-right changer					

EN 60601-1



ISO 13485:2016

ISO 14001:2015

OHSAS 18001:2007

ISO 9001:2015

MANUFACTURER COMPANY : BİLSER BİLİMSEL VE TEKNİK CİHAZLAR / ANKARA TURKEY

Tel:+90 312 395 06 23 *Fax: +90 312 395 06 73 *E-Mail: bilserbilsler@hotmail.com