

# Bilser®



*Since 1980, thousands of references ...*

## SURGICAL SUCTION UNIT



MODEL

BA-5



MODEL

BA-10

2 year guarantee

Bilser Brand Surgical Suctions, all transactions that may require surgical interventions and aspirator; operating room, intensive care, emergency units and is designed to be used in all clinics.

MODEL	Capacity	Jar	Jar Type	Pump	Vacuum	Type	Power
BA-5	5 lt	Single	Glass with safety float	Oil free	60 lt/min 600-650 mmHg  with vacuum adjustment control switch	Mobile  4 antistatic wheels, 2 with brakes	220-230 VAC 50/60 Hz
BA-10	2x5 lt	Double with left-right changer					

EN 60601-1



ISO 13485:2012

ISO 14001:2015

OHSAS 18001:2007

ISO 9001:2015

MANUFACTURER COMPANY : BİLSER BİLİMSEL VE TEKNİK CİHAZLAR / ANKARA TURKEY

Tel:+90 312 395 06 23 \*Fax: +90 312 395 06 73 \*E-Mail: bilserbilser@hotmail.com